



Elsaneek
Ref. **Absor-Protect 135**

QUIROTECH

**Absor
protect**
TECHNOLOGY

Verde/Green **3268C**

Azul/Blue **2727CP**

Azul/Blue **2717CP**



Composición / *Composition*

80% PES + 9% PA + 10% PU + 1% Antistatic (LATEX FREE)

Peso / *Weight* (Grams/m²)

135 ± 5%

Ligamento / *Ligament*

KNITED

UNE 40-339-76

Ancho / *Width* (cm)

150 ± 2

Encogimiento al lavado 75°C /
Washing shrinkage 75°C

3%

EN 6330/15797

Resist. a penetración de agua /
Water Penetration Resistance
(cm/H₂O)

> 1000

EN 20811

Absorción de agua / *Water absorbtion*
(mL/m²)

315

Regulación quirúrgica europea /
European surgical regulation

Pass
(55 washes and sterilizations)

EN13795

Antiestático / *Antistatic*

Pass

EN1149-3

Instrucciones de lavado / *Washing instructions:*



Lavar a 60°C - 75°C
(recomendado a 60°C)
Wash at 60°C - 75°C
(recommended at 60°C)



No permite el uso de lejía
Bleach NOT allowed



Secado: Secadora: de 60°C a 75°C máx. (según carga)
Esterilización: Vapor a 134°C / 4 minutos
Drying process: Tumble dryer: 60°C to 75°C max.
(depending on the load)
Sterilization: Vapour 134°C / 4 minutes



No planchar
Do not iron