

# MATTRESS AND PILLOW PROTECTOR

The most cost-effective, hygienic and safest solution on the market for hospitals and nursing homes



1

## Maximum absorption

High liquid absorption capacity (>1500 ml/m<sup>2</sup>) as a result of its high-capillarity technical fibres.



## Maximum waterproofness

Its Absorprotect technology guarantees maximum absorption and protection in accordance with the EN20811 standard.



## Image and comfort

Low-weight fabric for more convenient and quicker handling when making the bed, retains its whiteness and doesn't go yellow.

## Long-lasting

Suitable for intensive washes. Stays waterproof and retains absorption for up to 300 washes.



### Washing instructions:

Wash between 40°C and 80°C

Do not use bleach

Drying:  
Dryer: Max. temperature: 60°C - 80°C  
15-25 minutes (depending on load)

Do not iron